Precision balances KERN PBS · PBJ



Multifunctional laboratory balance with single-cell weighing system and EC type approval [M]



2 only PBS: **Adjusting program CAL** for quick setting of the balance accuracy, external test weights at an additional price



Simple recipeweighing and documenting with a combined tare/print function. In addition, the mixing ingredients for the recipe are numbered automatically and printed out with their corresponding number and weight value



Percentage determination: makes it possible to store a given weight value (100 %) and to determine deviations from this target value

Precision balances KERN PBS · PBJ







Features

- 11 only PBJ: Automatic internal adjustment in the case of a change in temperature and time-controlled at defined intervals
- · Metal housing: robust and sturdy
- Dosage aid: High-stability mode and other filter settings can be selected
- · Weighing with tolerance range (checkweighing): Input of an upper/lower limit value. A visual signal assists with portion division, dispensing or grading
- Totalizing of individual weighing results
- Identification number: 4 numeric positions, printed on calibration protocol
- Automatic data output to the PC/printer every time the balance is steady
- 3 Draft shield standard, on all models with readout [d] = 0,001 g, weighing space WxDxH 180x193x87 mm

Technical data

- Large backlit LCDdisplay, digit height 14 mm
- · Dimensions of weighing plate (stainless steel)
- A 108x105 mm
- B 170x180 mm, see larger picture
- Overall dimensions (without draft shield) WxDxH 209x322x78 mm
- Net weight approx. 3 kg
- Permissible ambient temperature 10 °C / 30 °C

Accessories

- Protective working cover over keyboard and housing, standard, can be reordered, for models with weighing plate sizes
- A KERN PBS-A01
- **B** KERN PBS-A02
- 4 Set for density determination for models with weighing plate sizes
- A KERN PBS-A04
- **B** KERN PBS-A03
- RS-232/Ethernet adapter to connect balances with an RS-232 interface to a network, using Ethernet, details see page 139, KERN YKI-01
- Suitable printers see page 138

- Single-cell advanced technology:
- Fully automatic manufactured weighing cell from one piece of material
- Stable temperature behavior
- Short stabilisation time: 3 sec
- Shock proof construction
- · High corner load performance

STANDARD



































| Model | Weighing | Readout | Verific. value | Minimum load | Repro- duci- | Linearity | Weighing plate | Options | | | |
|---------------|----------|---------|-------------------|-----------------|-----------------|-----------|----------------|--------------|--|-------------------------|---|
| | range | | | | | | | Verification | | DKD Calibr. Certificate | |
| | [Max] | [d] | [e] | [Min] | bility | | | M | | DKD | |
| KERN | g | g | g | g | g | g | | KERN | | KERN | |
| PBS 420-3M * | 420 | 0,001 | - | - | 0,001 | ± 0,002 | Α | - | | 963-127 | |
| PBS 620-3M | 620 | 0,001 | - | - | 0,001 | ± 0,002 | Α | - | | 963-127 | |
| PBS 4200-2M | 4200 | 0,01 | - | - | 0,01 | ± 0,02 | В | - | | 963-127 | |
| PBS 6200-2M | 6200 | 0,01 | - | - | 0,01 | ± 0,02 | В | - | | 963-128 | |
| PBS 8200-1M * | 8200 | 0,1 | - | - | 0,1 | ± 0,2 | В | - | | 963-128 | |
| | | • | | | | | | | | | • |
| PBJ 420-3M * | 420 | 0,001 | 0,01 | 0,02 | 0,001 | ± 0,002 | Α | 965-216 🗓 | | 963-127 | |
| PBJ 620-3M | 620 | 0,001 | 0,01 | 0,1 | 0,001 | ± 0,002 | Α | 965-201 🗓 | | 963-127 | |
| PBJ 4200-2M | 4200 | 0,01 | 0,1 | 0,5 | 0,01 | ± 0,02 | В | 965-216 🗓 | | 963-127 | |
| PBJ 6200-2M | 6200 | 0,01 | 0,1 | 1 | 0,01 | ± 0,02 | В | 965-202 🗓 | | 963-128 | |
| PBJ 8200-1M | 8200 | 0,1 | 1 | 5 | 0,1 | ± 0,2 | В | 965-217 🗓 | | 963-128 | |

Note: For applications that require verification, please order verification at the same time, initial verification at a later date is not possible. Verification at the factory, we need to know the full address of the location of use.

KERN Pictograms



Internal adjusting (CAL): Quick setting of the balance's accuracy with internal adjusting weight (motordriven).



Adjusting program (CAL): For quick setting of the balance's accuracy. External adjusting weight required.



Recipe level A: Separate memory for the weight of the tare container and the recipe ingredients (net total).



Rechargeable battery pack: rechargeable set.



RECIPE

Recipe level B: Internal memory for complete recipes with name and target value of the recipe ingredients. User guidance through displays.



Mains adapter: 230V/50Hz in standard version for Germany. On request GB, AUS or USA version.



Memory: Balance contains memories, e.g. for item data, weighing data, tare weights etc.



Data interface RS-232: To connect the balance to a printer, PC or network.



Recipe level C: Internal memory for complete recipes with name and target value of the recipe ingredients. User guidance through displays. Additional convenient functions, such as barcode and back calculation functions.



Power supply: integrated in balance. 230V/50Hz in Germany. More standards e. g. GB, AUS, USA on request.



Strain gauges: Electrical resistor on an elastic deforming body.



RS 232

RS 485 data interface: To connect the balance to a printer, PC or other peripheral devices. High tolerance against electromagnetic disturbance.



Percentage determination: Determining the deviation in % from the target value (100%).



Tuning fork principle: A resonating body is electromagnetically excited, causing it to oscillate.



USB data interface: To connect the balance to a printer, PC or other peripheral devices.



Weighing units: Can be switched to e. g. nonmetric units at the touch of a key. See balance model. Please refer to KERN's website for more details.



Electromagnetic force compensation: Coil in a permanent magnet. For the most accurate weighings.



Bluetooth data interface: To transfer data from the balance to a printer, PC or other peripheral devices.



Weighing with tolerance range: Upper and lower limiting can be programmed individually, e.g. dosing/sorting and portioning.



Single cell technology: Advanced version of the force compensation principle with the highest level of precision.



Control outputs (optocoupler, digital I/O): to connect relays, signal lamps, valves, etc.



Vibration-free weighing: (Animal weighing program) Vibrations are filtered out so that a stable weight is obtained.



Verification possible: The time required for verification is specified in the pictogram.



Interface for second balance: for direct connection of a second balance.



Spray and dust protection IPxx: The type of protection is shown by the pictogram. For details see the glossary.



DKD calibration possible: The time required for DKD calibration is shown in days in the pictogram.



Network interface: For connecting the scale to an Ethernet network. With KERN products you can also use a universal RS-232/LAN



Stainless steel: the balance is protected against corrosion.



Package shipment: The time required to manufacture the product internally is shown in days in the pictogram.



GLP/ISO record keeping: of weighing data with date, time and identification-no. Only with printers from KERN.



Suspended weighing: load support with hook on the underside of the balance.



Pallet shipment: The time required to manufacture the product internally is shown in days in the pictogram.



Piece counting: Reference quantities selectable. Display can be switched from piece to weight.



Battery operation: Ready for battery operation. The battery type is specified for each device.



Warranty: The warranty period is shown in the pictogram.

Precision is our business

To ensure the high precision of your balance KERN offers you the the appropriate test weight package for your balance, consisting of the test weight, box and DKD certificate, as proof of ist accuracy ... the best pre-requisite for proper balance

In the extensive KERN test weight range, you will find test weights in the international OIML error limit classes: E1, E2, F1, F2, M1, M2, M3 with weights from 1 mg - 2000 kg.

The KERN DKD calibration laboratory for electronic balances and weights has been accredited by DKD since 1994 and today is one of the most modern and best-equipped DKD calibration laboratories for balances, test weights and forcemeasurement in Europe.

(DKD = German Calibration Service)

Thanks to the high level of automation, we can carry out DKD calibration of balances, test weights and force-measuring devices 24 hours a day, 7 days a week.

Range of services:

- DKD calibration of balances with a maximum load of up to 6 t
- DKD calibration of weights in the range of 1 mg 500 kg
- · Database supported management of checking equipment and reminder service
- · Calibration of force-measuring devices
- DKD calibration certificates in the following languages D, GB, F, I, E, NL, PL

Do you have questions about your scale, the corresponsing test weight or the calibration service? Your KERN specialist dealer will be pleased to assist you.

Your KERN specialist dealer: